

SCHMITZ

A Sense of Confidence

varimed[®] **Furniture for OR theatres and outpatient departments**




varimed® OR theatre and outpatient department furniture



Model No.

Description


222.01.145



Surgeon's stool

5-arm pedestal in 18/10 stainless steel, antistatic plastic feet, rust-proof, seat disc in 18/10 stainless steel, height adjustment by means of threaded spindle with end stop, seat diam. 350 mm, height adjustment from 430 to 600 mm.


222.01.146



Surgeon's stool

as 222.01.145, but seat disc upholstered, antistatic cover. Height adjustment from 485 to 655 mm.


222.01.151



Surgeon's stool

5-arm pedestal in 18/10 stainless steel, antistatic double castors diam. 50 mm of plastic material, rust-proof, seat disc in 18/10 stainless steel, height adjustment by means of threaded spindle with end stop, seat diam. 350 mm, height adjustment from 430 to 600 mm.

222.01.153



Surgeon's stool


as 222.01.151, but seat disc upholstered, antistatic cover, black. Height adjustment from 485 to 655 mm.



Model No.

Description


222.01.165



Theatre seat

Specially-designed upholstered seat with antistatic cover, black, height adjustment from 560 to 790 mm by means of gas spring, with foot ring release. 5-arm pedestal powder coated in grey aluminium RAL 9007 with antistatic double plastic castors of 50 mm and circular black bumper.

222.01.167




Theatre seat

Specially-designed upholstered seat with antistatic cover, black, height adjustment from 560 to 790 mm by means of gas spring, with foot ring release, back rest with antistatic cover. 5-arm pedestal powder coated in grey aluminium RAL 9007 with antistatic double plastic castors of 50 mm and circular black bumper.




222.01.166



Theatre seat

Saddle-shaped pad with antistatic cover, black, height adjustment from 560 to 790 mm by means of gas spring, with foot ring release. 5-arm pedestal powder coated in grey aluminium RAL 9007 with antistatic double plastic castors of 50 mm and circular black bumper.

222.01.168



Theatre seat

as 222.166, but with back rest with antistatic cover, black. 5-arm pedestal powder coated in grey aluminium RAL 9007 with antistatic double plastic castors of 50 mm and circular black bumper.

varimed® OR theatre and outpatient department furniture




Please note:

The new version of instrument serving tables with hydraulic foot pump will be available as of January 2012. For orders placed before that date, the former version with ref. nos. 232.417 to 232.420 will still be supplied.

Model No.

Description

232.02.417



Instrument serving table C€/93/42/EEC

with foot-operated pump, height adjustment from 950 to 1,350 mm by means of foot bar, 18/10 stainless steel table top 700 x 500 mm, with raised edge, 360° rotatable, to be fixed in each position, rounded corners and edges; support column and frame in 18/10 stainless steel, pedestal clad in stainless steel, double plastic castors diam. 75 mm, ball bearing, antistatic, rust-proof.


Safe working load with table-top in any position: 30 kg.

232.02.418

Instrument serving table C€/93/42/EEC

as 232.02.417, but table top 600 x 400 mm.

232.02.419



Instrument serving table C€/93/42/EEC

as 232.02.417, but height adjustment from 800 to 1,200 mm

232.02.420


Instrument serving table C€/93/42/EEC

as 232.02.417, but table top 600 x 400 mm, height adjustment from 800 to 1,200 mm



232.02.417

232.02.426



Instrument serving table

18/10 stainless steel frame, 18/10 stainless steel table top 600 x 400 mm, non-rotatable, height adjustable from 900 to 1,320 mm, double plastic castors diam. 50 mm, antistatic, rust-proof. Safe working load: 12 kg.




232.02.426



Model No.

Description

232.02.427




Instrument serving table

manually operated height adjustment from 800 – 1,200 mm. Mobile 4-arm pedestal of 18/10 stainless steel and stainless steel table top with reflex free and disinfectant-proof surface. 4 swivel castors, resistant to abrasions, antistatic, diam. 75 mm, with friction bearing. Safe working load: 30 kg



232.02.427

232.02.446



Hand and arm operating table C€/93/42/EWG

pedestal of 18/10 stainless steel. Table top X-ray permeable, tiltable by means of ball-and-socket-joint, with clamping lever. Manually height adjustable from 650 – 950 mm. 4 swivel castors, lockable, resistant to abrasions, antistatic, diam. 75 mm, with friction bearing. Safe working load: 15 kg.

200.02.411

**Accessories:
Upholstered pad**

antistatic, X-ray permeable, detachable, 450 x 450 x 40 mm.




232.02.446
200.02.411

varimed[®] OR theatre and outpatient department furniture



Model No. Description

232.03.521




Instrument table

18/10 SST rectangular tube frame (frame also serves as protection rail), 18/10 SST table tops 600 x 450 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.522
232.03.523

Instrument table, as 232.03.521, but of 900 mm height
Instrument table, as 232.03.521, but of 1,000 mm height

232.03.531




Instrument table

18/10 SST rectangular tube frame, two stainless steel 18/10 shelves, top shelf 650 x 500 mm covering frame, lower one as inserted plate of 600 x 450 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.532
232.03.533

Instrument table, as 232.03.531, but of 900 mm height
Instrument table, as 232.03.531, but of 1,000 mm height

232.03.541




Instrument table

18/10 SST rectangular tube frame (frame also serves as protection rail), 18/10 SST table tops 700 x 600 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.542
232.03.543

Instrument table, as 232.03.541, but of 900 mm height
Instrument table, as 232.03.541, but of 1,000 mm height

232.03.551




Instrument table

18/10 SST rectangular tube frame, two stainless steel 18/10 shelves, top shelf 750 x 650 mm covering frame, lower one as inserted plate of 700 x 600 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.552
232.03.553

Instrument table, as 232.03.551, but of 900 mm height
Instrument table, as 232.03.551, but of 1,000 mm height

232.03.561




Instrument table

18/10 SST rectangular tube frame (frame also serves as protection rail), 18/10 SST table tops 900 x 600 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.562
232.03.563

Instrument table, as 232.03.561, but of 900 mm height
Instrument table, as 232.03.561, but of 1,000 mm height

232.03.571



Instrument table

18/10 SST rectangular tube frame, two stainless steel 18/10 shelves, top shelf 950 x 650 mm covering frame, lower one as inserted plate of 900 x 600 mm, detachable, table height 800 mm, 4 plastic double castors diam. 75 mm, antistatic, antirust, abrasion-proof, 2 of which lockable, with 4 bumpers. Safe working load per shelf: 50 kg

232.03.572
232.03.573

Instrument table, as 232.03.571, but of 900 mm height
Instrument table, as 232.03.571, but of 1,000 mm height

200.03.246

Accessories for instrument tables:

Drawer, all-stainless steel 18/10, effective size (w x d x h) 435 x 375 x 90 mm (not with 232.03.531 – 232.03.533, 232.03.551 – 232.03.553, 232.03.571 – 232.03.573 tables)

Swivel-mounted trays & holders, 1 round and 1 square stainless steel tray

Guard rail* 3-sided, stainless steel, for instrument tables 232.03.521 – 232.03.523

Guard rail* 3-sided, stainless steel, for instrument tables 232.03.541 – 232.03.543

Guard rail* 3-sided, stainless steel, for instrument tables 232.03.561 – 232.03.563


*subsequent attachment possible

200.03.132
200.03.185
200.03.186
200.03.187

Plastic castors with axle and fastening pin of stainless steel



242.04.125



Bowl stand

18/10 stainless steel frame, one 18/10 SST bowl, capacity 6 l, plastic double castors diam. 50 mm, antistatic, rust-proof, height 815 mm, **disassembled**.

Accessories:

Bowl, 18/10 stainless steel (extra)

200.04.307

242.04.135



Bowl stand


18/10 stainless steel frame, two 18/10 SST bowls, capacity 6 l ea., plastic double castors diam. 50 mm, antistatic, rust-proof, height 815 mm, **disassembled**.

Accessories:

Bowl, 18/10 stainless steel (extra)

200.04.307

242.04.131



Bowl stand CE/93/42/EEC

for heating of physiological saline solutions, adjustment knob to switch the temperature from 35° C to max. 60° C, with safety thermostat and scale reading. Teflon-coated heating plate for easy cleaning. 5-arm base of stainless steel, height adjustable by servo-assisted gas cylinders; plastic castors diam. 75 mm, rust-proof. 2 castors with brake, 3 antistatic, SST bowl with litre-scale, capacity 5 l. Earth contact. Voltage: 230 V, 50 Hz, 400 W, connection cable 5 m with plug.

Accessories:

5 L Bowl, stainless steel 18/10, with litre-scale (extra)
Fastening ring, for aseptic cover

200.04.560
200.04.561

Important: For reasons of safety only use original bowls from Schmitz.


varimed® OR theatre and outpatient department furniture



Model No.

Description


242.05.235



Drip-feed stand CE/93/42/EEC

18/10 stainless steel frame, disassembled, 5-arm base, upper section with 4 hooks, drip glass, double plastic castors diam. 50 mm, antistatic, rust-proof; height adjustment from 1,300 to 2,120 mm, one-hand action. Safe working load: 2 kg (1 hook each).

251.05.101



Step

18/10 stainless steel frame, fine-ribbed anti-slip tread in 18/10 stainless steel, floor pads antistatic, tread height 240 mm.


Step, as 251.05.101, but tread height 120 mm.

Step, as 251.05.101, but tread height 300 mm.

Safe working load: 185 kg

251.05.103
251.05.104

251.05.102




Step

18/10 stainless steel frame, fine-ribbed anti-slip treads in 18/10 stainless steel, floor pads antistatic.

Safe working load: 185 kg

252.05.118



Refuse bucket

removable 12 l bucket in 18/10 stainless steel, 18/10 stainless steel frame with profiled plastic wrap-around bumper, double plastic castors diam. 50 mm, antistatic, rust-proof; frame height 230 mm.

252.05.119

Refuse bucket

as 252.05.118 above, with pushing handle in 18/10 stainless steel.

200.05.137


200.05.135

Accessories:

Bucket lid in 18/10 stainless steel

12 L Bucket in 18/10 stainless steel (extra)

252.05.130



Soiled linen container

18/10 stainless steel frame, 1 laundry bag holder, firmly welded, with rubber ring, 4 swivel castors, diam. 75 mm.

200.05.231

200.05.232

200.05.176


Accessories:

Lid plastic, with fastening attachment, w/ plastic handle

Lid plastic, without fastening attachment, w/ plastic handle

Linen bag self-emptying

252.05.135



Soiled linen container

18/10 stainless steel frame, 2 laundry bag holders, firmly welded, with rubber ring, 4 swivel castors, diam. 75 mm.

200.05.231

200.05.232

200.05.176

Accessories:

Lid plastic, with fastening attachment, with plastic handle

Lid plastic, without fastening attachment, with plastic handle
















Linen bag self-emptying

varimed®

OR theatre and outpatient department furniture

Accessories

Material description

Model No.	Description	Model No.	Description	Model No.	Description
200.03.130	Round SST tray with holder of SST, for instrument tables	200.03.131	Rectangular SST tray with holder of SST, for instrument tables	200.03.132	Swivel-mounted trays and holders 1 round and 1 rectangular steel tray, for instrument tables
					
200.03.301	Tray holder round, 18/10 SST, for 200.03.303	200.03.302	Tray holder rectangular, 18/10 SST, for 200.03.304	200.03.303	SST tray round
					
200.03.304	SST tray rectangular	200.05.135	12 L Bucket SST 18/10 for 252.05.118 – 252.05.119	200.05.137	Lid for 12 l. bucket, 18/10 SST
					
200.05.176	Linen bag for 252.05.130 and 252.05.135	200.05.231	Lid plastics, with fixation and plastic handle, for 252.05.130 and 252.05.135	200.05.232	Lid plastics, w/o fixation, plastic handle, for 252.05.130 and 252.05.135
					
200.03.185	Guard rail 18/10 SST, for instrument tables 232.03.521 – 232.03.523	200.03.186	Guard rail 18/10 SST, for instrument tables 232.03.541 – 232.03.543	200.03.187	Guard rail 18/10 SST, for instrument tables 232.03.561 – 232.03.563
					
200.03.246	Drawer completely made of 18/10 stainless steel, effective size 435 x 375 x 90 mm (w x d x h) for instrument tables	200.04.307	Bowl for bowl stands 242.04.125 and 242.04.135, 18/10 SST, diam. of 320 mm, depth of 130 mm	200.04.560	Bowl stainless steel 18/10, with litre-scale for 242.04.131
					
200.04.561	Fastening ring for aseptic cover for 242.04.131	200.02.411	Pad antistatic, X-ray permeable, for hand and arm operating table 232.02.446		
					

Material description

Frames in SST, round or rectangular tubing, Material 4301, electropolished. Stainless steel tops as well as claddings are of ground, brushed finish (grain size 180). All castors are resistant to abrasions. Seat upholsteries of top quality, hard-wearing foam rubber with antistatic black cover.



Schmitz u. Söhne GmbH & Co. KG

P. O. Box 14 61
58734 Wickede (Ruhr)
Germany (postal address)
Zum Ostenfeld 29
58739 Wickede (Ruhr)
Germany (visiting address)
Phone +49 (0) 2377-84-0
Fax +49 (0) 2377-84-135
www.schmitz-soehne.com
export@schmitz-soehne.de

Technical service, hotline:
Phone +49 (0) 2377 84549
Fax +49 (0) 2377 84210
service@schmitz-soehne.de

All articles including corresponding accessories marked with the CE label meet the standards for medical devices acc. to MPG and the EC regulation 93/42/EEC.



The products may only be used for their intended purpose of use. Please also observe our operating manual!

Schmitz abroad: Switzerland

SCHMITZ AG
Schloss-Strasse 67
8207 Schaffhausen
Phone +41 (0) 52 6 43 25 05
Fax +41 (0) 52 6 43 26 91
f.poncet@bluewin.ch

SCHMITZ AG
François Poncet
Chemin de Bois 32
1020 Renens VD
Phone +41 (0) 21 6 91 74 66
Fax +41 (0) 21 6 91 74 62
Mobile 079-2128334
f.poncet@bluewin.ch

Spain

SCHMITZ u. Söhne Ibérica S.L.
Centro de Negocios Tartessos
Edificio Artemisa
C/Pollensa, n° 2 – Oficina 3
28290 Las Rozas (Madrid)
Phone +34 (0) 91 6 40 45 14
Fax +34 (0) 91 6 36 63 05
schmitz@schmitz.es

Italy

SCHMITZ ITALIA S.r.l.
Daniela Capeletti
Piazza Garibaldi 20
15076 Ovada (AL)
Phone +39 (01 43) 8 03 50
Fax +39 (01 43) 83 42 53
Mobile 349-6338723
info@schmitz-italia.it
www.schmitz-italia.it

France

SCHMITZ FRANCE SAS
27, rue du Lieutenant Colonel Prévost
69006 Lyon
Phone +33 (0) 3 83 22 94 67
Fax +33 (0) 3 83 22 87 66
info@schmitz-france.fr

Russia and Eastern Europe

Schmitz u. Söhne GmbH & Co. KG
Sergey Kosolapov
Bornaische Str. 180
04279 Leipzig
Germany
Phone +49 (0) 3 41 3 37 80 34
Fax +49 (0) 3 41 3 37 80 35
Mobile +49 (0) 1 52 0 90 77 05 8
kosolapov@schmitz-soehne.de

China

Schmitz u. Söhne GmbH & Co. KG
Guangzhou Representative Office
World Trade Centre,
North Tower 3105
Huan Shi East Road 371-375
Guangzhou 510095
Phone +86 (0) 20 87 61 33 01
Fax +86 (0) 20 87 60 39 23
info@schmitz-soehne.cn

Latin America

Schmitz u. Söhne GmbH & Co. KG
Alfredo Gentile
Av. Córdoba 1513 Piso 8, Oficina A
Buenos Aires, Argentina
Phone +54 (0)11 4811 4642
Fax +54 (0)11 4811 4642
Mobile 0911 6887 2086
gentile@schmitz-soehne.de



The complete Schmitz-products range



Gynaecology,
Urology and
Proctology
The perfect
equipment to meet
all requirements



DIAMOND
A class of its own



OPX Operating
Tables for Out-
patient Departments
and Operating
Theatres
Unique Flexibility
and Reliability



STL 285, STS 282,
STX 280
The Patient Transport
Systems for Emer-
gencies, Outpatient
Depts., Intensive
Care and X-Ray



varimed®
The trolleys –
our experts for
everyday
clinical practice



varimed®
Medical Furniture
for Examination and
Treatment Areas



Partura®
The Schmitz
delivery table
comfortable and
safe delivery –
with focus on pre-
ventive obstetrics