

Transfer mattress for easy and comfortable transfer of patients for patient stretchers STS 282 and STL 285



80 mm thick, of resistant, hard-wearing PU material, antistatic, black, suitable for X-ray imaging, underside of high quality special sliding material, 1,950 x 600 mm, with 4 looped handles along each side. To be secured on both sides of the stretcher by Velcro straps.



Transfer mattress, Model No. 206.389
for patient stretcher STS 282.101

Transfer mattress, Model No. 207.389
for patient stretcher STL 285.080 and 285.100

Getting to the point – Your key benefits



Secured easily and safely by Velcro straps on both sides



- High sliding quality and optimum comfort
- Less strenuous and easy to use
- Helps to save time and costs as it can be handled by just 1 person
- Ideal for transfer to X-ray tables
- Suitable for patients up to 225 kg
- No additional investment, as it can also be used as a standard mattress

Important:

The transfer mattress may only be used on Schmitz patient stretchers 282.101, 285.080 and 285.100 manufactured after 02/2009!

The transfer mattress cannot be secured on previous models nor is it suitable for use on patient stretchers from other manufacturers!