

varimed® Height adjustable Examination and Massage couches




NEW



- Large and ergonomic height adjustment range
- No displacement of table-top during height adjustment
- Moulded base provides generous legspace




varimed Examination and Massage couch CE



Examination and Massage couch Model no. 211.410

hydraulic height adjustment 500 to 920 mm by means of foot lever. Servo-assisted head section adjustment from -20° to $+40^{\circ}$. 2-section upholstery 70mm thick, leatherette cover acc. to Schmitz colour chart. Powder coated chassis in white aluminium RAL 9006 or grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover. With 4 levelling feet, not electrically conductive. Safe working load: 225 kg
Available from May 2011



Examination and Massage couch Model no. 211.420

as 211.410, but electromotive height adjustment from 500 to 920 mm by means of foot control unit, 2.40 m connecting cable with mains plug. Upon request, hand control instead of foot control unit.
Available now

Accessories:

Mobility Model no. 200.025

4 closed double castors \varnothing 100 mm, 2 of which el. conductive, central lock per axis, 1 castor with directional lock

Safety ON/OFF switch Model no. 200.030

turns the foot control unit on or off.

Note!

the safety ON/OFF switch 200.030 is required by the law in Germany, Austria and Switzerland.

Paper roll holder Model no. 200.093

for paper roll up to 590 mm wide, mounting at foot end only

Paper roll Model no. 200.103

500 mm wide

Head section with nose outlet Model no. 200.050

Leatherette comfort

	mint green 80
	moss green 81
	caribbean blue 82
	atlantic blue 83
	blossom white 84
	parchment beige 85
	gold yellow 86
	tobacco brown 87
	granitgrau 89
	granite grey 90
	titan grey 91
	red carmine 92
	natural brown 93

Leatherette classic

	white 53
	silver grey 55
	ocean 52
	black 57